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
The following pressing buttons have a double function depending how they are used or needed:

If pressed once: **goes up**
 If kept pressed: **It allows you to set program 1**

If pressed once: **Goes down**
 If kept pressed: **It allows you to set program 2**

If pressed once: **It Works as "Enter" key.**
 If kept pressed: **It allows you to set program 3**

“FlexTone” Thermo Injection Programming or reprogramming:

Step 1: Plug equipment in to 110 voltz as well as air intake to 115 psi (8bar).

Step 2: Press any of the following buttons (being careful machine is off), equipment will make a noise:

PROGRAM 1 / PROGRAM 2 / PROGRAM 3

Step 3: Release the button. You will observe in the screen the following data:

Program:1

Temp:_____ Time:_____ Cooling:_____

On the left side of screen an arrow will show on, with said arrow select the value to be modified with buttons “down” or “up”.

Step 4: Once you are positioned on the value to be modified, press the “Enter” button. Arrow will change its shape.

Program 1:

Temp:_____ Time:_____ Cooling:_____

Using buttons “down” or “up” you will be able to increase or decrease the value according to your requirements.

Step 5: Press “Enter” in order to save the new value.

Step 6: If you require it, please repeat this operation with each of the values of the program.

Step 7: Once program changes have been performed, turn the equipment off. In this way you will be saving the new value configuration. You will be able to verify new programming by turning the machine on. If required, make same operation with programs 2 and 3.

Please take in consideration that you can reprogram your equipment anytime you wish. It is recommended having a program with minimum values in order to make the piston head change each time it is needed. Head change is made by using 26g. cartridges in order to avoid material waste).

How to use Flex Tone thermo injection equipment.

Step 1: Plug equipment to 110 voltz and also air intake to 115 psi (8 bar).

Step 2: By using the “power” button which is located on the left side, turn your equipment on.

Step 3: Select your program 1, 2 or 3 by pressing “up” and “down” buttons.

Step 4: Once you have selected your program, press “Enter”, a red light will turn on, it indicates selected program has been started. Resistance begins to heat at that moment.

Step 5: Wait until program reaches and stabilizes its temperature. A Green light will show meaning system is ready to receive another cartridge.

Step 6: Lubricate cartridge to be used and insert it at the oven. press the “ok” button. Equipment will start working on programmed temperature and time.

Step 7: As soon as machine has reached programmed temperature, piston will automatically be activated changing its light to blue. Piston will keep pressure during programmed cooling time.

Step 8: Once cooling time has elapsed, piston will automatically retract, until then, you will be able to take the muffle out from the equipment. Light will change to green color.

Step 9: In order to clean or expulse injected cartridge debris left inside the oven, pressure “ok” button without releasing it, then pressure “expulse” button. Once this operation has been done, equipment will be ready to start a new cycle. Light will keep Green.

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Maintenance:

In order to maintain your equipment in good conditions for a long working period, it is recomendable to clean plunger each time it has been used. This way. With this, you will avoid metal and material debris accumulation on each injection. Provide a wider cleaning to plunger each 30 days preferably .

Steps in order to provide an adequate Maintenance:

- Repeat from step 3
- Only for Maintenance Purposes, we will provide minimum values as temperature, time and cooling that are instructed for equipment.
- When equipment is in cooling process, we will disable air supply, the piston will keep elevated.
- Turn the equipment off.
- Once machine is cool, we will take out piston head by using the right tooling.
- Once it has been taken, we will be able to give a better cleaning within that area as well as to plunger by helping us with Flex Tone release spray.
- Once it is clean, place head on its original position and secure it with proper tools.
- Connect air supply, piston will retract to its original position.
- Reprogram if necessary.
- You can start working with your equipment normally.
- By following this steps we can also exchange different head caliber, depending on the size of cartridge to be used.

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